

Product Information

MemDX™ Membrane Protein Human TMEM198 (Transmembrane protein 198) Full Length

Cat. No.: **MPC1442K**

This product is for research use only and is not intended for diagnostic use.

This product is a 39.4 kDa Human TMEM198 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM198

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

39.4 kDa

TMD

7

Sequence

MPGTVATLRFQLLPPEPDDAFWGAPCEQPLERRYQALPALVCIMCCLFGV
VYCFFGYRCFAVLFLTGLLFGSVVIFLLCYRERVLETQLSAGASAGIAL
GIGLLCGLVAMLVRSVGLFLVGLLLGLLAAAAALLGSAPYYQPGSVWGPL
GLLLGGLLCALLTLRWPRPLTTLATAVTGAALIATAADYFAELLLLGRY
VVERLRAAPVPPLCWRSWALLALWPLLSLMGVLVQWRVTAEGDSHTEVVI
SRQRRRVQLMRIRQQEDRKEKRRKKRPPRAPLRGPRAPPRPGPPDPAYRR
RPVPIKRFNGDVLSPSYIQSFRDRQTGSSLSSFMASPTDADYEYGSRGPL
TACSGPPVRV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM198

Full Name

Transmembrane protein 198

Alternative Names

TMEM198; TMEM198A; Transmembrane protein 198

Gene ID

[130612](#)

UniProt ID

[Q66K66](#)