

Product Information

MemDX™ Membrane Protein Human TMEM179 (Transmembrane protein 179) Full Length

Cat. No.: **MPC1444K**

This product is for research use only and is not intended for diagnostic use.

This product is a 26.3 kDa Human TMEM179 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM179

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

26.3 kDa

TMD

4

Sequence

MALNNFLFAQCACYFLAFLFSFVVVPLSENGHDFRGRCLLFTEGMWLSA
NLTVQERERFTVQEHWGPPAACRFSLASLLSLLAAAHAWRTLFFLCKGH
EGSFFSAFLNLLVSAFVFLVFIASITIVSVGFTMWCDTITEKGTVPHSCE
ELQDIDLELGVDNSAFYDQFAIAQFGLWASWLAWLAITTLAFLKVYHNYR
QEDLLDSLIHEKELLLARPSRPSFQEEKSAVI

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM179

Full Name

Transmembrane protein 179

Alternative Names

TMEM179; C14orf90; TMEM179A; transmembrane protein 179A; Transmembrane protein 179

Gene ID

[388021](#)

UniProt ID

[Q6ZVK1](#)