

Product Information

MemDX™ Membrane Protein Human VSTM4 (V-set and transmembrane domain containing

4) Full Length

Cat. No.: **MPC1460K**

This product is for research use only and is not intended for diagnostic use.

This product is a 36.1 kDa Human VSTM4 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

VSTM4

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

36.1 kDa

TMD

1

Sequence

MRLLALAAAALLARAPAPEVCAALNVTVSPGPVVDYLEGENATLLCHVSQ
KRRKDSLLAVRWFFAHSFDSQEALMVKMTKLRVVQYYGNFSRSKRRLR
LLEEQRGALYRLSVLTLQPSDQGHYVCRVQEISRHRNKWTAWSNGSSATE
MRVISLKASEESSFEKTKETWAFFEDLYVYAVLVCCVGILSILLFMLVIV
WQSVFNKRKSRVRHYLVKCPQNSSGETVTSVTS LAPLQPKKGKRQKEKPD
IPPAVPAKAPIPTFHKPKLLKQQRKVTLPKIAEENLTYAELELIKPHRA
AKGAPTSTVYAQILFEENKL

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

VSTM4

Full Name

V-set and transmembrane domain containing 4

Alternative Names

VSTM4; C10orf72; V-set and transmembrane domain-containing protein 4; V-set and transmembrane domain containing 4

Gene ID

[196740](#)

UniProt ID

[Q8IW00](#)