

Product Information

MemDX™ Membrane Protein Human TMEM251 (Transmembrane protein 251) Full Length

Cat. No.: **MPC1461K**

This product is for research use only and is not intended for diagnostic use.

This product is a 18.7 kDa Human TMEM251 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM251

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

18.7 kDa

TMD

2

Sequence

MPKPPDYSELSDSLTLAVGTGRFSGPLHRAWRMMNFRQRMGWIGVGLYLL
ASAAAFYYVFEISETYNRLALEHIQQHPPEPLEGTTWTHSLKAQLLSLPF
WVWTVIFLVPYLQMFLFLYSCTRADPKTVGYCIIPICLAVICNRHQAFVK
ASNQISRLQLIDT

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

TMEM251

Full Name

Transmembrane protein 251

Alternative Names

TMEM251; C14orf109; UPF0694 transmembrane protein C14orf109; Transmembrane protein 251

Gene ID

[26175](#)

UniProt ID

[Q8N6I4](#)