

Product Information

MemDX™ Membrane Protein Human LAX1 (Lymphocyte transmembrane adaptor 1) Full

Length

Cat. No.: **MPC1462K**

This product is for research use only and is not intended for diagnostic use.

This product is a 44 kDa Human LAX1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

LAX1

Protein Length

Full length

Protein Class

Immunity

Molecular Weight

44 kDa

TMD

1

Sequence

MDGVTPTLSTIRGRTLESSTLHVTPRSLDRNKDQITNIFSGFAGLLAILL
VVAVFCILWNWNKRKKRQVPYLRVTVMPLLTLPQTRQRAKNIYDILPWRQ
EDLGRHESRSMRIFSTESLLSRNSESPEHVPSQAGNAFQEHTAHIHATEY
AVGIYDNAMVPQMCGNLTPSAHCINVRASRDCASISSEDSHDYVNVPTAE
EIAETLASTKSPSRNLFVLPSTQKLEFTTEERDEGCGDAGDCTSLYSPGAE
DSDSLSNREGSSQISNDYVNMTGLDLSAIQERQLWVAFQCCRDIYENVPAA
DPSGSQQQAEKDVPSSNIGHVEDKTDDPGTHVQCVRKRTFLASGDYADFQP
FTQSEDSQMKHREEMSNESSDYENVLTAKLGGRDSEQGPQTQLLPDE

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

LAX1

Full Name

Lymphocyte transmembrane adaptor 1

Alternative Names

LAX1; LAX; LAT-like membrane associated protein; linker for activation of x cells; membrane-associated adapter protein
LAX; Lymphocyte transmembrane adaptor 1

Gene ID

[54900](#)

UniProt ID

[Q8IWW1](#)