

Product Information

MemDX™ Membrane Protein Human TMEM200A (Transmembrane protein 200A) Full Length

Cat. No.: **MPC1482K**

This product is for research use only and is not intended for diagnostic use.

This product is a 54.3 kDa Human TMEM200A membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM200A

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

54.3 kDa

TMD

2

Sequence

MIATGGVITGLAALKRQDSARSQQHVNLSPPATQEKKPIRRRPRADV
VRGKIRLYSPSGFFLILGVLISIIGIAMAVLGYWPQKEHFIDAETTLSTN
ETQVIRNEGGVVVRRFFEQLHSDKMKMLGPFTMGIGIFIFICANAILHEN
RDKETKIIHMRDIYSTVIDIHTLRIKEQRQMNGMYTGLMGETEVKQNGSS
CASRLAANTIASFSGFRSSFRMDSSVEEDELMLNEGKSSGHLMPPLSDS
SVSVFGLYPPPSKTTDDKTSGSKKCETKSIVSSSISAFSLPVIKLNNCVI
DEPSIDNITEDADNLKSRSRNLSMDSLVVPLPNTSESFQPVSTVLPNNNS
IGESLSSQYKSSMALGPGAGQLLSPGAARRQFGSNTSLHLLSSHSKSLDL
DRGPSTLTVQAEQRKHPSWPRLDRNNSKGYMKLENKEDPMDRLLVPQVAI
KKDFTNKEKLLMISRSHNNLSFEHDEFSLNNLKRGTSERF

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM200A

Full Name

Transmembrane protein 200A

Alternative Names

TMEM200A; TTMA; TTMC; KIAA1913; two transmembrane C; HBE61; Transmembrane protein 200A

Gene ID

[114801](#)

UniProt ID

[Q86VY9](#)