

Product Information

MemDX™ Membrane Protein Human TM4SF18 (Transmembrane 4 L six family member 18)

Full Length

Cat. No.: **MPC1484K**

This product is for research use only and is not intended for diagnostic use.

This product is a 22.2 kDa Human TM4SF18 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TM4SF18

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

22.2 kDa

TMD

4

Sequence

MGSRKCGGCLSCLLIPLALWSIIVNILLYFPNGQTSYASSNKLNTNYVWYF
EGICFSGIMMLIVTTVLLVLENNNNYKCCQSENC SKKYVTLLSIIFSSLG
IAFSGYCLVISALGLVQGPYCRTL DGWEYAFEGTAGRFLTDSSIWIQCLE
PAHVVEWNIILFSILITLSGLQVIICLIRVVMQLSKILCGSYSVIFQPGI
I

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TM4SF18

Full Name

Transmembrane 4 L six family member 18

Alternative Names

TM4SF18; L6D; transmembrane 4 L6 family member 18; Transmembrane 4 L six family member 18

Gene ID

[116441](#)

UniProt ID

[Q96CE8](#)