

# Product Information

## MemDX™ Membrane Protein Human TM4SF18 (Transmembrane 4 L six family member 18)

### Full Length

Cat. No.: **MPC1484K**

This product is for research use only and is not intended for diagnostic use.

This product is a 22.2 kDa Human TM4SF18 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

TM4SF18

#### Protein Length

Full length

#### Protein Class

Transporter

#### Molecular Weight

22.2 kDa

#### TMD

4

#### Sequence

MGSRKCGGCLSCLLIPLALWSIIVNILLYFPNGQTSYASSNKLNTNYVWYF  
EGICFSGIMMLIVTTVLLVLENNNNYKCCQSENC SKKYVTLLSIIFSSLG  
IAFSGYCLVISALGLVQGPYCRTL DGWEYAFEGTAGRFLTDSSIWIQCLE  
PAHVVEWNIILFSILITLSGLQVIICLIRVVMQLSKILCGSYSVIFQPGI  
I

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements

**Form**

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

TM4SF18

**Full Name**

Transmembrane 4 L six family member 18

**Alternative Names**

TM4SF18; L6D; transmembrane 4 L6 family member 18; Transmembrane 4 L six family member 18

**Gene ID**

[116441](#)

**UniProt ID**

[Q96CE8](#)