

Product Information

MemDX™ Membrane Protein Human TMEM14A (Transmembrane protein 14A) Full Length

Cat. No.: **MPC1511K**

This product is for research use only and is not intended for diagnostic use.

This product is a 10.7 kDa Human TMEM14A membrane protein expressed in *E.coli*-based cell-free expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM14A

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

10.7 kDa

TMD

3

Sequence

MDLIGFGYAALVTFGSIFGYKRRGGVPSLIAGLFVGCLAGYGAYRVSNDK
RDVKVSLFTAFFLATIMGVRFKRSKKIMPAGLVAGLSLMMILRLVLLLL

Product Description

Expression Systems

E.coli-based cell-free expression system

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM14A

Full Name

Transmembrane protein 14A

Alternative Names

TMEM14A; PTD011; C6orf73; Transmembrane protein 14A

Gene ID

[28978](#)

UniProt ID

[Q9Y6G1](#)