

# Product Information

## MemDX™ Membrane Protein Human TMEM65 (Transmembrane protein 65) Full Length

Cat. No.: **MPC1530K**

This product is for research use only and is not intended for diagnostic use.

This product is a 25.4 kDa Human TMEM65 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

TMEM65

#### Protein Length

Full length

#### Protein Class

Transporter

#### Molecular Weight

25.4 kDa

#### TMD

3

#### Sequence

MSRLLPLLRSRTARSLRPGPAAAAAPRPPSWCCCGRLLALAPPGGLPGG  
PRRLGTHPKKEPMEALNTAQGARDFIYSLHSTERSCLLKELHRFESIAIA  
QEKLEAPPPTPGQLRYVFIHNAIPFIGFGFLDNAMIVAGTHIEMSIGII  
LGISTMAAAALGNLVSDLAGLGLAGYVEALASRLGLSIPDLTPKQVDMWQ  
TRLSTHLGKAVGVTIGCILGMFPLIFFGGGEEDEKLETKS

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements

**Form**

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

TMEM65

**Full Name**

Transmembrane protein 65

**Alternative Names**

TMEM65; Transmembrane protein 65

**Gene ID**

[157378](#)

**UniProt ID**

[Q6PI78](#)