

Product Information

MemDX™ Membrane Protein Human TMEM129 (Transmembrane protein 129, E3 ubiquitin ligase) Full Length

Cat. No.: **MPC1553K**

This product is for research use only and is not intended for diagnostic use.

This product is a 40.4 kDa Human TMEM129 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM129

Protein Length

Full length

Protein Class

Transferase

Molecular Weight

40.4 kDa

TMD

3

Sequence

MDSPEVTFTLAYLVFAVCFVFTPNEFHAAGLTVQNLLSGWLGSDDAAFPV
FHLRRTAATLLCHSLLPLGYYVGMCLAASEKRLHALSQAPEAWRLFLLLA
VTLPSIACILIIYYWSRDRWACHPLARTLALYALPQSGWQAVASSVNTFR
RIDKFATGAPGARVIVTDTWVMKVTTYRVHVAQQQDVHLTVTESRQHELS
PDSNLPVQLLTIRVASTNPAVQAFDIWLNSTEYGELCEKLRAPIRRAHV
VIHQSLGDLFLETFAVLVEVNPAYSVPSSELEACIGCMQTRASVKLVKT
CQEAATGECQQCYCRPMWCLTCMGKWFASRQDPLRPDTWLASRVPCPTCR
ARFCILDVCTVR

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM129

Full Name

Transmembrane protein 129, E3 ubiquitin ligase

Alternative Names

TMEM129; D4S2561E; E3 ubiquitin-protein ligase TM129; RING-type E3 ubiquitin transferase TM129; transmembrane protein 129, E3 ubiquitin protein ligase; Transmembrane protein 129, E3 ubiquitin ligase

Gene ID

[92305](#)

UniProt ID

[A0AVI4](#)