

Product Information

MemDX™ Membrane Protein Human TMEM116 (Transmembrane protein 116) Full Length

Cat. No.: **MPC1574K**

This product is for research use only and is not intended for diagnostic use.

This product is a 27.4 kDa Human TMEM116 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM116

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

27.4 kDa

TMD

4

Sequence

MKHTQSGQSTSPLVIDYTCRVQMAFVFSSLIPLLLMTPVFCLGNTSECF
QNFSQSHKCILMHSPPSAMAELPPSANTSV CSTLYFYGIAIFLGSFVLSL
LTIMVLLIRAQTLTKPKQDTKLHMAFYVLQALTATSQGLLNCGVYGWTQHKFH
GPAVILMIKLTQPKQDTKLHMAFYVLQALTATSQGLLNCGVYGWTQHKFH
QLKQEARRDADTQTPLLC SQKRFYSRGLNSLESTLTFPASTSTIF

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM116

Full Name

Transmembrane protein 116

Alternative Names

TMEM116; Transmembrane protein 116

Gene ID

[89894](#)

UniProt ID

[Q8NCL8](#)