

# Product Information

## MemDX™ Membrane Protein Human TMEM95 (Transmembrane protein 95) Full Length

Cat. No.: **MPC1591K**

This product is for research use only and is not intended for diagnostic use.

This product is a 19.5 kDa Human TMEM95 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

TMEM95

#### Protein Length

Full length

#### Protein Class

Transporter

#### Molecular Weight

19.5 kDa

#### TMD

1

#### Sequence

MWRLALGGVFLAAAQACVFCRLPAHDLSGRLARLCSQMEARQKECGASPD  
FSAFALDEVSMNKVTEKTHRVLRVMEIKEAVSSLPSYWSWLRKTKLPEYT  
REALCPPACRGSTTLNCSTCKGTEVSCWPRKRCFPGSQDLWEAKILLLS  
IFGAFLLLGVLSLLVESHHLQAKSGL

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements

#### Form

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

TMEM95

**Full Name**

Transmembrane protein 95

**Alternative Names**

TMEM95; UNQ9390; sperm-egg fusion protein TMEM95; Transmembrane protein 95

**Gene ID**

[339168](#)

**UniProt ID**

[Q3KNT9](#)