

Product Information

MemDX™ Membrane Protein Human TMEM185B (Transmembrane protein 185B) Full Length

Cat. No.: **MPC1609K**

This product is for research use only and is not intended for diagnostic use.

This product is a 40.5 kDa Human TMEM185B membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM185B

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

40.5 kDa

TMD

7

Sequence

MNPRGLFQDFNPSKFLIYTCLLLFSVLLPLRLDGIIQWSYWAVFAPIWLW
KLLVVAGASVGAGVWARNPRYRTEGEACVEFKAMLIAVGIHLLLMFEVL
VCDRVERGTHFWLLVFMPPLFFVSPVSAACVWGFRHDSLELEILCSVNI
LQFIFIALKLDRIIHWPLVVFVPLWILMSFLCLVLYYIVWSLLFLRSL
DVVAEQRRTHVTMAISWITIVVPLLTFEVLLVHRLDGHNTFSYVSIFVPL
WLSLLTLMATTFRRKGGNHWWFGIRRDQCFLLEIFPFLREYGNISYDLH
HEDSEDAEETSVPEAPKIAPIFGKKARVVITQSPGKYVPPPKLNIDMPD

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM185B

Full Name

Transmembrane protein 185B

Alternative Names

TMEM185B; FAM11B; ee3_2; family with sequence similarity 11, member B; putative transmembrane protein 185B; transmembrane protein 185B (pseudogene); Transmembrane protein 185B

Gene ID

[79134](#)

UniProt ID

[Q9H7F4](#)