

Product Information

MemDX™ Membrane Protein Human TMEM156 (Transmembrane protein 156) Full Length

Cat. No.: **MPC1614K**

This product is for research use only and is not intended for diagnostic use.

This product is a 34.3 kDa Human TMEM156 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM156

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

34.3 kDa

TMD

2

Sequence

MTKTALLKLFVAIVITFILPEYFKTPKERTLELSCLEVCLQSNFTYSL
SSLNFSFVTLQPVRETQIIMRIFLNPSNFRNFTRTCQDITGEFKMCSSC
LVCESKGNMDFISQEQTSKVLIRRGSM EVKANDFHSPCQHFNFSVAPLVD
HLEEYNTTCHLKNHTGRSTIMEDEPSKEKSINYTCRIMEYPNDCIHISLH
LEMDIKNITCSMKITWYILVLLVFIFLIILTIRKILEGQRRVQKWQSHRD
KPTSVLLRGSDSEKLRLALNVQVLSAETTQRLPLDQVQEVLPPIPEL

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM156

Full Name

Transmembrane protein 156

Alternative Names

TMEM156; Transmembrane protein 156

Gene ID

[80008](#)

UniProt ID

[Q8N614](#)