

Product Information

MemDX™ Membrane Protein Human VSTM1 (V-set and transmembrane domain containing

1) Full Length

Cat. No.: **MPC1688K**

This product is for research use only and is not intended for diagnostic use.

This product is a 26.1 kDa Human VSTM1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

VSTM1

Protein Length

Full length

Protein Class

Immunity

Molecular Weight

26.1 kDa

TMD

1

Sequence

MTAEFLSLLCLGLCLGYEDEKKNEKPPKPSLHAWPSSVVEAESNVTLCQ
AHSQNVTFVLRKVNDSGYKQEQSSAENAEFPFTDLKPKDAGRYFCAYKT
TASHEWSESSEHLQLVVTDKHDELEAPSMKTDTRTIFVAIFSCISILLF
LSVFIYRCSQHSSSSEESTKRTSHSKLPEQEAADLSNMERVSLSTAD
PQGVTYAELSTSALSEAASDTTQEPPGSHEYAALKV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

VSTM1

Full Name

V-set and transmembrane domain containing 1

Alternative Names

VSTM1; SIRL1; SIRL-1; UNQ3033; V-set and transmembrane domain-containing protein 1; LAIR homolog; OSCAR-like transcript-1; signal inhibitory receptor on leukocytes-1; V-set and transmembrane domain containing 1

Gene ID

[284415](#)

UniProt ID

[Q6UX27](#)