

Product Information

MemDX™ Membrane Protein Human TEDDM1 (Transmembrane epididymal protein 1) Full Length

Cat. No.: **MPC1739K**

This product is for research use only and is not intended for diagnostic use.

This product is a 31.3 kDa Human TEDDM1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TEDDM1

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

31.3 kDa

TMD

6

Sequence

MILKGCLLYPLCSPRNKQRCARLWKIAYGGLLKIVTGSLLTFYVVLCLDG
GMVLMRKQVPSRFMYPKEWQHLLTMFILLTLNGCVDFMSKNVLPQRCVGLE
KGTLLVLIIEYELLLMVSHVKDSEGVELHVYSLILVVFLLLLVLTAEWA
PNMCHLQLMETFLILMMGSWLMQAGFILYRPVSGYPWQDDDISDIMFVTT
FFCWHVMINASFLGIIYGFSFVYHCFRPSLKLTPKEAPYYASTPGPLY
KLLQEVEQSEKEDQALLLPKSSP

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TEDDM1

Full Name

Transmembrane epididymal protein 1

Alternative Names

TEDDM1; Epididymal protein 9; epididymis secretory sperm binding protein Li 45e; human epididymis-specific protein 9; putative membrane protein HE9; transmembrane protein 45C; Transmembrane epididymal protein 1

Gene ID

[127670](#)

UniProt ID

[Q5T9Z0](#)