

Product Information

MemDX™ Membrane Protein Human LDAF1 (Lipid droplet assembly factor 1) Full Length

Cat. No.: **MPC1752K**

This product is for research use only and is not intended for diagnostic use.

This product is a 17.5 kDa Human LDAF1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

LDAF1

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

17.5 kDa

TMD

4

Sequence

MAKEEPQSI~~SRDLQELQKKLSLLIDSFQNN~~SKV~~VAFMKSPVGQYLD~~SHPF
LAFTLLVFIVMSAVPVGFFLLIV~~LTTLAALLGVIILEGLVISVGGF~~SLL
CILCGLGFVSLAMSGMMIASYVVVSSLISCWFSPRPLTQQNTSCDFLPAM
KSAEFEGLYQE

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

LDAF1

Full Name

Lipid droplet assembly factor 1

Alternative Names

LDAF1; TMEM159; PROMETHIN; transmembrane protein 159; Lipid droplet assembly factor 1

Gene ID

[57146](#)

UniProt ID

[Q96B96](#)