

Product Information

MemDX™ Membrane Protein Human TMEM88B (Transmembrane protein 88B) Full Length

Cat. No.: **MPC1754K**

This product is for research use only and is not intended for diagnostic use.

This product is a 17 kDa Human TMEM88B membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM88B

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

17 kDa

TMD

2

Sequence

MSEQGRETEEEEGGGGASDTAPMLPRGPPDHQASALTCPGWSGPPLLPGR
LLAGLLHLLLPAAAFLLVLLPAAAVVYLGFLCHSRVHPAPGPRCRAF
DRGSAALIVFGLLSLPPLLVLASAVRARLARLRLLPPPAGTPGPRRPP
GRPDEDEQLCAWV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM88B

Full Name

Transmembrane protein 88B

Alternative Names

TMEM88B; Transmembrane protein 88B

Gene ID

[643965](#)

UniProt ID

[A6NKF7](#)