

Product Information

MemDX™ Membrane Protein Human SPPL3 (Signal peptide peptidase like 3) Full Length

Cat. No.: **MPC1774K**

This product is for research use only and is not intended for diagnostic use.

This product is a 42.2 kDa Human SPPL3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SPPL3

Protein Length

Full length

Protein Class

Protease

Molecular Weight

42.2 kDa

TMD

9

Sequence

MAEQTYSWAYSLVDSSQVSTFLISILLIVYGSFRSLNMDFENQDKEKDSN
SSSGSFNGNSTNNSIQTIDSTQALFLPIGASVLLVMFFFFDSVQVVFTI
CTAVLATIAFAFLLPMCQYLTRPCSPQNKISFGCCGRFTAELLFSLS
VMLVLIWVLTGHWLLMDALAMGLCVAMIAFVRLPSLKVSCLLLSGLLIYD
VFWVFFSAYIFNSVMVKVATQPADNPLDVLNRKHLGPNVGRDVPRLSL
PGKLVFPSSTGSHFSMLGIGDIVMPGLLLCFVLRDNYKKQASGDSCGAP
GPANISGRMQKVSYPFHCTLIGYFVGLLTATVASRIHRAAQPALLYLVPFT
LLPLLTMAYLKGDLLRMWSEPFHKSSSSRFLEV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SPPL3

Full Name

Signal peptide peptidase like 3

Alternative Names

SPPL3; IMP2; PSH1; PSL4; PRO4332; MDHV1887; signal peptide peptidase-like 3; SPP-like 3; intramembrane protease 2; presenilin homologous protein 1; presenilin-like protein 4; Signal peptide peptidase like 3

Gene ID

[121665](#)

UniProt ID

[Q8TCT6](#)