

Product Information

MemDX™ Membrane Protein Human TMEM249 (Transmembrane protein 249) Full Length

Cat. No.: **MPC1798K**

This product is for research use only and is not intended for diagnostic use.

This product is a 27 kDa Human TMEM249 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM249

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

27 kDa

TMD

2

Sequence

MPKGRAGSLPTTSIGWRFQLWFLGLTCPERHLARRLKNNSFYPPFVQQEPN
VVFLEYLDLTLWKGMLLFIISVVLVSFSSLREVQKQETWVFLVYGVGVL
WLVISSLPRRRLVLNHTRGVYHFSIQGRTVCQGPHLVYVRLALSSDAHG
RCFFHLVLGGHRLEPLVLVQLSEHYEQMEYLGRIYARKLNINYFDYLATS
YRHVVRHWPPPGAGTVMGKSPMGHKPSSSQSSLEV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM249

Full Name

Transmembrane protein 249

Alternative Names

TMEM249; C8ORFK29; putative transmembrane protein C8orfK29; Transmembrane protein 249

Gene ID

[340393](#)

UniProt ID

[Q2WGJ8](#)