

# Product Information

## MemDX™ Membrane Protein Human TMEM14EP (Transmembrane protein 14E, pseudogene)

### Full Length

Cat. No.: **MPC1967K**

This product is for research use only and is not intended for diagnostic use.

This product is a 13.6 kDa Human TMEM14EP membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

TMEM14EP

#### Protein Length

Full length

#### Protein Class

Transporter

#### Molecular Weight

13.6 kDa

#### TMD

2

#### Sequence

MQMDPGPQVPLYWLGFVYAALAALGGISGYAKVGSVQSPSAGFFFSELAG  
LDASQPSRNPKEHLSSPVYIWDLARYYANKILTLWNIYACGFSCRCLLIV  
SKLGSMYGEQILSVVAMSQGLGLMKN

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements

**Form**

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

TMEM14EP

**Full Name**

Transmembrane protein 14E, pseudogene

**Alternative Names**

TMEM14EP; TMEM14E; Transmembrane protein 14E, pseudogene

**Gene ID**

[645843](#)

**UniProt ID**

[Q6UXP3](#)