

Product Information

MemDX™ Membrane Protein Human RHBDL3 (Rhomboid like 3) Full Length

Cat. No.: **MPC1971K**

This product is for research use only and is not intended for diagnostic use.

This product is a 45.2 kDa Human RHBDL3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

RHBDL3

Protein Length

Full length

Protein Class

Protease

Molecular Weight

45.2 kDa

TMD

7

Sequence

MGEHPSPGPAVAACAEAEERIEELEPEAEERLPAAPEDHWKVLFDQFDPGN
TGYISTGKFRSLLESHSSKLDPHKREVLLALADSHADGQIGYQDFVSLMS
NKRSNSFRQAILQGNRRLSSKALLEEKGLSLSQRLIRHVAYETLPREIDR
KWYYDSYTCPPPWFMITVTLLLEVAFFLYNGVSLGQFVLQVTHPRYLKNS
LVYHPQLRAQVWRYLTYIFMHAGIEHLGLNVVLQLLVGVPLEMVHGATRI
GLVYVAGVVAGSLAVSVADMTAPVVGSSGGVYALVSAHLANIVMNWSGMK
CQFKLLRMAVALICMSMEFGRAVWLRFHPSAYPPCPHPSFVAHLGGVAVG
ITLGVVVLRNIEQRLQDQSLWWIFVAMYTVFVLFVAVFWNIFAYTLLDLKL
PPPP

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

RHBDL3

Full Name

Rhomboid like 3

Alternative Names

RHBDL3; VRHO; RHBDL4; rhomboid-related protein 3; rhomboid, veinlet-like 3; rhomboid, veinlet-like 4; ventrhoid transmembrane protein; Rhomboid like 3

Gene ID

[162494](#)

UniProt ID

[P58872](#)