

Product Information

MemDX™ Membrane Protein Human TLCD4 (TLC domain containing 4) Full Length

Cat. No.: **MPC2001K**

This product is for research use only and is not intended for diagnostic use.

This product is a 30 kDa Human TLCD4 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TLCD4

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

30 kDa

TMD

6

Sequence

MEINTKLLISVTCISFFTFQLLFYFVSYWFSKAVSPGFNSLSFKKKIEWN
SRVVSTCHSLVVGIFGLYIFLFDEATKADPLWGGPSLANVNIAIASGYLI
SDLSIIILYWKVIGDKFFIMHHCASLYAYYLVKNGVLAYIGNFRLLAEL
SSPFVNQRWFFEALKYPKFSKAIVINGILMTVVFVIRIASMLPHYGFMY
SVYGTEPYIRLGVLQLSWVISCVLDVMNVMWMIKISKGCIVISHIRQ
EKAKNSLQNGKLD

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TLCD4

Full Name

TLC domain containing 4

Alternative Names

TLCD4; TMEM56; TLC domain-containing protein 4; transmembrane protein 56; TLC domain containing 4

Gene ID

[148534](#)

UniProt ID

[Q96MV1](#)