

Product Information

MemDX™ Membrane Protein Human TSPAN33 (Tetraspanin 33) Full Length

Cat. No.: **MPC2044K**

This product is for research use only and is not intended for diagnostic use.

This product is a 31.5 kDa Human TSPAN33 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TSPAN33

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

31.5 kDa

TMD

4

Sequence

MARRPRAPAASGEEFSFVSPLVKYLLFFFNMLFWVISMVMVAVGVYARLM
KHAEAAALACLAVDPAILLIVGVLMFLLTFCGCIGSLRENICLLQTFSLC
LTAVFLLQLAAGILGFVFSKARGKVSEIINNAIVHYRDDLDLQNLIDFG
QKKFSCCGGISYKDWSONMYFNCSNPSRERCSVPYSCCLPTPDQAVIN
TMCQGGMQAFDYLEASKVIYTNCGIDKLVNWIHSNLFLLGGVALGLAIPQ
LVGILLSQILVNQIKDQIKLQLYNQQHRADPWY

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TSPAN33

Full Name

Tetraspanin 33

Alternative Names

TSPAN33; PEN; PEN.; TSPAN-33; tetraspanin-33; penumbra; proerythroblast new membrane; Tetraspanin 33

Gene ID

[340348](#)

UniProt ID

[Q86UF1](#)