

Product Information

MemDX™ Membrane Protein Human EMC6 (ER membrane protein complex subunit 6) Full Length

Cat. No.: **MPC2079K**

This product is for research use only and is not intended for diagnostic use.

This product is a 12 kDa Human EMC6 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

EMC6

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

12 kDa

TMD

3

Sequence

MAAVVAKREGPPFISEAAVRGNAAVLDYCRTSVSALSGATAGILGLTGLY
GFIFYLLASVLLSLLLILKAGRRWNKYFKSRRPLFTGGLIGGLFTYVLFW
TFLYGMVHVY

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -70°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

EMC6

Full Name

ER membrane protein complex subunit 6

Alternative Names

EMC6; TMEM93; RAB5IFL; transmembrane protein 93; ER membrane protein complex subunit 6

Gene ID

[83460](#)

UniProt ID

[Q9BV81](#)