

Product Information

MemDX™ Membrane Protein Human TRGC2 (T cell receptor gamma constant 2) Full Length

Cat. No.: **MPC2136K**

This product is for research use only and is not intended for diagnostic use.

This product is a 21.6 kDa Human TRGC2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TRGC2

Protein Length

Full length

Protein Class

Receptor; Immunity

Molecular Weight

21.6 kDa

TMD

1

Sequence

DKQLDADVSPKPTIFLPSIAETKLQKAGTYLCLLEKFFPDIIKIHWQEKK
SNTILGSQEGNTMKTNDTYMKFSWLTVP EESLDKEHRCIVRHENKNGID
QEIIFFPIKTDVTTVDPKYNYSKDANDVITMDPKDNWSKDANDTLLLQLT
NTSAYTYLLLLLKS VVYFAITCCLLRRTAFCCNGEKS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TRGC2

Full Name

T cell receptor gamma constant 2

Alternative Names

TRGC2; TCRGC2; TRGC2(2X); TRGC2(3X); T-cell receptor, gamma, constant region C2; T cell receptor gamma constant 2

Gene ID

[6967](#)

UniProt ID

[P03986](#)