

Product Information

MemDX™ Membrane Protein Human TMDD1 (Transmembrane and death domain 1) Full

Length

Cat. No.: **MPC2149K**

This product is for research use only and is not intended for diagnostic use.

This product is a 33.5 kDa Human TMDD1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMDD1

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

33.5 kDa

TMD

1

Sequence

MAARTLASALVLTWVWALAPAGAVDAMGPHAAVRLAELLTPEECGHFRS
LLEAPEPDVEAELSRLSEDRLARPEPLNTTSGSPSRRRRREAAEDPAGRV
AGPGEVSDGCREALAAWLAPQAASLSWDRLARALRRSGRPDVARELGKNL
HQQATLQLRKFGQRFLPRPGAAARVPFAPAPRPRRAAVPAPDWDALQLIV
ERLPQPLYERSPMGWAGPLALGLLTGFGALGTGALVVLLTLWITGGDGD
RASPGSPGPLATVQGWETKLLLPKERRAPPGAWAADGPDSPSPHSALAL
SCKMGAQSWGSGALDGL

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMDD1

Full Name

Transmembrane and death domain 1

Alternative Names

TMDD1; C12orf81; hCG1648055; transmembrane and death domain protein 1; Uncharacterized protein C12orf81; Transmembrane and death domain 1

Gene ID

[112163659](#)

UniProt ID

[P0DPE3](#)