

Product Information

MemDX™ Membrane Protein Human COX8C (Cytochrome c oxidase subunit 8C) Full Length

Cat. No.: **MPC2160K**

This product is for research use only and is not intended for diagnostic use.

This product is a 8.1 kDa Human COX8C membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

COX8C

Protein Length

Full length

Protein Class

Oxidoreductase

Molecular Weight

8.1 kDa

TMD

1

Sequence

MPLLRGRCPARRHYRRLALLGLQPAPRFAHSGPPRQRPLSAAEMAVGLV
FFTTFLTPAAYVLGNLKQFRRN

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

COX8C

Full Name

Cytochrome c oxidase subunit 8C

Alternative Names

COX8C; COX8-3; cytochrome c oxidase subunit 8C, mitochondrial; COX VIII-3; cytochrome c oxidase polypeptide 8; cytochrome c oxidase polypeptide VIII; cytochrome c oxidase subunit 8-3; cytochrome c oxidase subunit VIII; cytochrome c oxidase subunit VIIIC; Cytochrome c oxidase subunit 8C

Gene ID

[341947](#)

UniProt ID

[Q7Z4L0](#)