

Product Information

MemDX™ Membrane Protein Human STMP1 (Short transmembrane mitochondrial protein 1)

Full Length

Cat. No.: **MPC2197K**

This product is for research use only and is not intended for diagnostic use.

This product is a 5.2 kDa Human STMP1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

STMP1

Protein Length

Full length

Protein Class

Immunity

Molecular Weight

5.2 kDa

TMD

1

Sequence

MLQFLLGFTLGNNVGMVLAQNYDIPNLAKKLEEEIKKDLDAKKKPPSA

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

STMP1

Full Name

Short transmembrane mitochondrial protein 1

Alternative Names

STMP1; Mm47; C7orf73; PL-5283; uncharacterized protein C7orf73; Short transmembrane mitochondrial protein 1

Gene ID

[647087](#)

UniProt ID

[E0CX11](#)