

Product Information

MemDX™ Membrane Protein Human TOMM20 (Translocase of outer mitochondrial membrane 20) Full Length

Cat. No.: **MPC2203K**

This product is for research use only and is not intended for diagnostic use.

This product is a 16.2 kDa Human TOMM20 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TOMM20

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

16.2 kDa

TMD

1

Sequence

MVGRNSAIAAGVCGALFIGYCIYFDRKRRSDPNFKNRLRERRKKQKLAKE
RAGLSKLPDLKDAEAVQKFFLEEIQLGEELLAQGEYEKGV DHLTNAI AVC
GQPQQLLQVLQQTLP PPVFQMLLTKLPTISQRIVSAQSLAEDDVE

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TOMM20

Full Name

Translocase of outer mitochondrial membrane 20

Alternative Names

TOMM20; MAS20; MOM19; TOM20; mitochondrial import receptor subunit TOM20 homolog; mitochondrial 20 kDa outer membrane protein; outer mitochondrial membrane receptor Tom20; translocase of outer mitochondrial membrane 20 homolog type II; Translocase of outer mitochondrial membrane 20

Gene ID

[9804](#)

UniProt ID

[Q15388](#)