

Product Information

MemDX™ Membrane Protein Human MXRA8 (Matrix remodeling associated 8) Full Length

Cat. No.: **MPC2215K**

This product is for research use only and is not intended for diagnostic use.

This product is a 49.1 kDa Human MXRA8 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

MXRA8

Protein Length

Full length

Protein Class

Cell adhesion

Molecular Weight

49.1 kDa

TMD

1

Sequence

MALPSRILLWKLVLQSSAVLLHSGSSVPAAAGSSVSESAVSWEAGARA
VLRCQSPRMVWTQDRLHQRVLHWDLRGPGGGPARRLLDLYSAGEQRVY
EARDRGRLELSASAFDDGNFSLIRAVEETDAGLYTCNLHHHYCHLYESL
AVRLEVTDGPPATPAYWDGEKEVLAVARGAPALLTCVNRGHVWTDHVEE
AQQVHWDQRPPGVPHDRADRLDLYASGERRAYGPLFLRDRVAVGADAF
ERGDFSLRIEPLVADEGTYSCHLHHHYCGLHERRVFHLTVAEPHAEP
RGSPGNGSSSHSGAPGPDPTLARGHNVINVIVPESRAHFFQQLGYVLATLL
LFILLVTVLLAARRRRGGYEYSDQKSGKSKGKDVNLAEFAVAAGDQMLY
RSEDIQLDYKNNILKERAELAHSPAPAKYIDLDKGFRKENCK

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

MXRA8

Full Name

Matrix remodeling associated 8

Alternative Names

MXRA8; ASP3; matrix remodeling-associated protein 8; limitrin; matrix-remodeling-associated protein 8; matrix-remodelling associated 8; Matrix remodeling associated 8

Gene ID

[54587](#)

UniProt ID

[Q9BRK3](#)