

Product Information

MemDX™ Membrane Protein Human TMEM191C (Transmembrane protein 191C) Full Length

Cat. No.: **MPC2229K**

This product is for research use only and is not intended for diagnostic use.

This product is a 39.2 kDa Human TMEM191C membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM191C

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

39.2 kDa

TMD

1

Sequence

MCRATLGLPLPPIVIQPARRSLPPIVTPASRRLGPRGGRHLGSVSTAMAA
TQELLLQLQKDNRDGRQRKQELEKLMRGLEAESESLNQRLQDLSEERERSL
LRRRSQAAQPLQGEAREAAERERAERVRRRLEEAERHKEYLEQHSRQLQEQ
WEELSSQLFYYGGELQSQKSTEQQLAAQLVTLQNELELAETKCALQEEKL
QQDALQTAEAWAIFQEQTVVVLQVRPHSDAKVPPASPPDLGRCDGQLRGV
QYSTESLMEEMARADRETRLFGGPRALAIRRCVLGALQVLLTLPLFLGL
SLLWTVLLDPGAVSAWLWSLTSETTLRRLRYTSLPLELRANGLLPT

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TMEM191C

Full Name

Transmembrane protein 191C

Alternative Names

TMEM191C; KB-1323B2.6; Transmembrane protein 191C

Gene ID

[645426](#)

UniProt ID

[A6NGB0](#)