

# Product Information

## MemDX™ Membrane Protein Human TMEM275 (Transmembrane protein 275) Full Length

Cat. No.: **MPC2234K**

This product is for research use only and is not intended for diagnostic use.

This product is a 17.2 kDa Human TMEM275 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

TMEM275

#### Protein Length

Full length

#### Protein Class

Transporter

#### Molecular Weight

17.2 kDa

#### TMD

2

#### Sequence

MPPAEKSEGPPVPAPAERARGRVPGLPSPALCCACGLCALLAGVNVTLAG  
AFASFLPEHNALLVGLALLVLALGFFAACCVCSRRLAPRGRSAAAAGP  
GQGGGRAGPVALEMESSEPTAQDTTAVQLSPAVSAASSGCSSPGPSPLAL  
EAPAPAAVLCALRSEGVQLNPPRARAAP

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements

#### Form

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

TMEM275

**Full Name**

Transmembrane protein 275

**Alternative Names**

TMEM275; Transmembrane protein 275

**Gene ID**

[105378696](#)

**UniProt ID**

[A0A0U1RQS6](#)