

Product Information

MemDX™ Membrane Protein Human OR4D1 (Olfactory receptor family 4 subfamily D member 1) Full Length

Cat. No.: **MPC2246K**

This product is for research use only and is not intended for diagnostic use.

This product is a 35.2 kDa Human OR4D1 membrane protein expressed in Baculovirus/Insect expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

OR4D1

Protein Length

Full length

Protein Class

GPCR

Molecular Weight

35.2 kDa

TMD

7

Sequence

MEPQNTTQVSMFVLLGFSQTQELQKFLFLLVYVTTIVGNLLIMVTVT
FDCRLHTPMYFLLRNALIDLCYSTVTSPKMLVDFLHETKTISYQGCMQAQ
IFFFHLLGGGTVFVFLSVMAYDRYIAISQPLRYVTIMNTQLCVGLVVAWV
GGFVHSIVQLALILPLPFCGPNILDNFYCDVPQVLR LACTDTSLLEFLMI
SNSGLLVIIWFLLLISYTVILVMLRSHSGKARRKAASTCTTHIIVSMI
FIPCIYITWPFPTFLMDKAVSISYVTMPMLNPMIYTLRNQDMKAAMRR
LGKCLVICRE

Product Description

Expression Systems

Baculovirus/Insect expression system

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

OR4D1

Full Name

Olfactory receptor family 4 subfamily D member 1

Introduction

Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms.

Alternative Names

OR4D1; OR4D3; OR4D4P; TPCR16; OR17-23; olfactory receptor 4D1; olfactory receptor 4D3; olfactory receptor OR17-23 pseudogene; olfactory receptor TPCR16; olfactory receptor, family 4, subfamily D, member 3; olfactory receptor, family 4, subfamily D, member 4 pseudogene; seven transmembrane helix receptor; Olfactory receptor family 4 subfamily D member 1

Gene ID

[26689](#)

UniProt ID

[Q15615](#)