

Product Information

MemDX™ Membrane Protein Human PVRIG (PVR related immunoglobulin domain containing) Full Length

Cat. No.: **MPC2247K**

This product is for research use only and is not intended for diagnostic use.

This product is a 34.3 kDa Human PVRIG membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

PVRIG

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

34.3 kDa

TMD

3

Sequence

MRTEAQVPALQPPEPGLEGAMGHRTLVLPPWLLTLCVTAGTPEVWVQVRM
EATELSSFTIRCGFLGSGSISLVTWSWGGPNGAGGTTLAVLHPERGIRQW
APARQARWETQSSISLILEGSGASSPCANTTFCKFASFPEGSWEACGSL
PPSSDPGLSAPPTPAPILRADLAGILGVSGVLLFGCVYLLHLLRRHKHRP
APRLQPSRTSPQAPRARAWAPSQASQAALHVPYATINTSCRPATLDTAHP
HGGPSWWASLPTHAAHRPQGPAAWASTPIPARGSFVSVENGLYAQAGERP
PHTGPGLTLFPDPRGPRAMEGPLGVR

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

PVRIG

Full Name

PVR related immunoglobulin domain containing

Alternative Names

PVRIG; CD112R; C7orf15; transmembrane protein PVRIG; CD112 receptor; nectin-2 receptor; poliovirus receptor related immunoglobulin domain containing; poliovirus receptor-related immunoglobulin domain-containing protein; PVR related immunoglobulin domain containing

Gene ID

[79037](#)

UniProt ID

[Q6DKI7](#)