

Product Information

MemDX™ Membrane Protein Human OST4 (Oligosaccharyltransferase complex subunit 4, non-catalytic) Full Length

Cat. No.: **MPC2251K**

This product is for research use only and is not intended for diagnostic use.

This product is a 4.1 kDa Human OST4 membrane protein expressed in *E.coli*. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

OST4

Protein Length

Full length

Protein Class

Transporter

Molecular Weight

4.1 kDa

TMD

1

Sequence

MITDVQLAIFANMLGVSLFLLVVLYHYVAVNNPKKQE

Product Description

Expression Systems

E.coli

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

OST4

Full Name

Oligosaccharyltransferase complex subunit 4, non-catalytic

Alternative Names

OST4; dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit 4; oligosaccharyltransferase 4 homolog; Oligosaccharyltransferase complex subunit 4, non-catalytic

Gene ID

[100128731](#)

UniProt ID

[P0C6T2](#)