

Product Information

MemDX™ Membrane Protein SARS2 S Mutation (S477R) (Surface glycoprotein) Full Length

Cat. No.: **MPC2304K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order SARS2 S Mutation (S477R) membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

SARS2

Target Protein

S Mutation (S477R)

Protein Length

Full length

Protein Class

Host-virus interaction

TMD

1

Sequence

MFVFLVLLPLVSSQCVNLTTTRTQLPPAYTNSFTRGVYYPDKVFRSSVLHS
TQDLFLPFFSNVTWFHAIHVSGTNGTKRFDNPVLPFNDGVYFASTEKSN
IRGWIFGTTLDSTQSLIVNNATNVVIVKVFCEQFCNDPFLGVYYHKNNK
SWMESEFRVYSSANNCTFEYVSQPFLMDLEGKQGNFKNLREFVFKNIDGY
FKIYSKHTPINLVRDLPGGFSALEPLVDLPIGINITRFQTLALHRSYLT
PGDSSSGWTAGAAAYVGYLQPRTFLLKYNENGTITDAVDCALDPLSETK
CTLKSFTVEKGIYQTSNFRVQPTESIVRFPNITNLCPFGEVFNATRFASV
YAWNKRISNVCVADYSVLVNSASFSTFKCYGVSPTKLNDLCFTNVYADSF
VIRGDEVQRQIAPGQTGKIADYNYKLPDDFTGCVIAWNSNNLDSKVGGNYN
YLRLFRKSNLKPFERDISTEIQAGRTPCNGVEGFNCYFPLQSYGFQPT
NGVGYQPYRVVLSFELLHAPATVCGPKKSTNLVKNKCVNFNFNGLTGTG
VLTESNKKFLPFQQFGRDIADTTDAVRDPQTLEILDITPCSFGGVSVITP
GTNTSNQVAVLYQDVNCTEVPVAIHADQLTPTWRVYSTGSNVFQTRAGCL
IGAETHVNNSECDIPIGAGICASYQTQTNSPRRARSVASQSIIAYTMSLG
AENSVAYSNNISIAIPTNFTISVTTEILPVSMTKTSVDCTMYICGDSTEC
NLLLQYGSFCTQLNRALTGIAVEQDKNTQEVFAQVKQIYKTPPIKDFGGF
NFSQILPDPSKPSKRSFIEDLLFNKVTADAGFIKQYGDCLGDIAARDLI
CAQKFNGLTVLPLLTDEMIQAQYTSALLAGTITSGWTFGAGAALQIPFAM
QMAYRFNGIGVTVQNVLYENQKLIANQFNSAIGKIQDSLSTASALGKLQD
VVNQNAQALNTLVKQLSSNFGAIVSVLNDILSRDLKVEAEVQIDRLITGR
LQSLQTYVTQQLIRAAEIRASANLAATKMSECVLGQSKRVDFCGKGYHLM
SFPQSAPHGVVFLHVTVYVPAQEKNFTTAPAICHGDKAHFPREGVFVSNGT
HWFVTQRNFYEPQIITTDNTFTVSGNCDVVIGIVNNTVYDPLQPELDSFKE
ELDKYFKNHTSPDVLGDISGINASVVNIQKEIDRLNEVAKNLNESLIDL

QELGKYEQYIKWPWYIWLGFIAGLIAIVMTIMLCCMTSCCCLKGCCSC
GSCCKFDEDDSEPVKGVKLHYT

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

S Mutation (S477R)

Full Name

Surface glycoprotein

Introduction

Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) is an enveloped, positive-sense, single-stranded RNA virus that causes coronavirus disease 2019 (COVID-19). Virus particles include the RNA genetic material and structural proteins needed for invasion of host cells. Once inside the cell the infecting RNA is used to encode structural proteins that make up virus particles, nonstructural proteins that direct virus assembly, transcription, replication and host control and accessory proteins whose function has not been determined.~ The structural proteins of SARS-CoV-2 include the envelope protein (E), spike or surface glycoprotein (S), membrane protein (M) and the nucleocapsid protein (N). The spike glycoprotein is found on the outside of the virus particle and gives coronavirus viruses their crown-like appearance. This glycoprotein mediates attachment of the virus particle and entry into the host cell. S protein is an important target for vaccine development, antibody therapies and diagnostic antigen-based tests.

Alternative Names

S Mutation (S477R); structural protein; spike protein; Surface glycoprotein

Gene ID

[43740568](#)

UniProt ID

[P0DTC2](#)