

Product Information

MemDX™ Membrane Protein Human FCRL6 (Fc receptor like 6) Full Length

Cat. No.: **MPC2431K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human FCRL6 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

FCRL6

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MLLWTAVLLFVPCVGKTVWLYLQAWPNPVFEGDALTLRCQGWKNTPLSQV
KFYRDGKFLHFSKENQTL SMGAATVQSRGQYSCSGQVMYIPQTFTQTSET
AMVQVQELFPPPVLSAIPSPPEPREGSLVTLRCQTKLHPLRSALRLLFSFH
KDGH TLQDRGPHPELCIPGAKEGDSGLYWCEVAPEGGQVQKQSPQLEVRV
QAPVSRPVLTLHHGPADPAVGDMVQLLCEAQRGSPPILYSFYLDEKIVGN
HSAPCGGTTSLLPVKSEQDAGNYSCEAENSVSRRERSEPKKLSLKGSQVL
FTPASNWLVPWLPASLLGLMVIAAALLVYVRSWRKAGPLPSQIPPTAPGG
EQCPLYANVHHQKGKDEGVVYSVVRHRTSKRSEARSAEFTVGRKDSSIICA
EVRCLQPSEVSSTEVNMR SRLQEPLSDCEEVLC

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

FCRL6

Full Name

Fc receptor like 6

Alternative Names

FCRL6; FcRH6; Fc receptor-like protein 6; Fc receptor-like protein 7; IFGP6; IgSF type I transmembrane receptor; fc receptor homolog 6; fcR-like protein 6; leukocyte receptor; Fc receptor like 6

Gene ID

[343413](#)

UniProt ID

[Q6DN72](#)