

Product Information

MemDX™ Membrane Protein Human C16orf54 (Chromosome 16 open reading frame 54) Full Length

Cat. No.: **MPC2444K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human C16orf54 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

C16orf54

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MPLTPEPPSGRVEGPPAWEAAPWPSLPCGPCIPIMLVLATLAALFILTTA
VLAERLFRRALRPDP SHRAPTLVWRPGGELWIEPMGTARERSEDWYGSAV
PLLTDRAPEPPTQVGTLEARATAPPAPSAPNSAPSNLGPQTVLEVPARST
FWGPQPWEGRPATGLVSWAEPEQRPEASVQFGSPQARRQRP GSPDPEWG
LQPRVTLEQISAFWKREGRTSVGF

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

C16orf54

Full Name

Chromosome 16 open reading frame 54

Alternative Names

C16orf54; transmembrane protein C16orf54; Chromosome 16 open reading frame 54

Gene ID

[283897](#)

UniProt ID

[Q6UWD8](#)