

Product Information

MemDX™ Membrane Protein Human TIMM23B (Translocase of inner mitochondrial membrane 23 homolog B) Full Length

Cat. No.: **MPC2456K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human TIMM23B membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TIMM23B

Protein Length

Full length

Protein Class

Transporter

TMD

2

Sequence

MEGGGGSGDKTTGVLGFFGAGEAGYSHADLAGVPLTG MNPLCPYLNVDP
RYLVQDTDEFILPTGANKTWGRFELAFFTIGGCCMTGA AFGAMNGLRLGL
KETQNM AWSKPGNVQILNMVTRQGALWANTLGSLALLYSAFGVIEKTRG
AEDDLNTVAAGTMTGMLYKCTVSEMALDSPFCVLLSGS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TIMM23B

Full Name

Translocase of inner mitochondrial membrane 23 homolog B

Alternative Names

TIMM23B; mitochondrial import inner membrane translocase subunit Tim23B; putative mitochondrial import inner membrane translocase subunit Tim23B; Translocase of inner mitochondrial membrane 23 homolog B

Gene ID

[100652748](#)

UniProt ID

[Q5SRD1](#)