

Product Information

MemDX™ Membrane Protein Human MPZL1 (Myelin protein zero like 1) Full Length

Cat. No.: **MPC2467K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human MPZL1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

MPZL1

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MAASAGAGAVIAAPDSRRWLWSVLAAALGLLTAGVSALEVYTPKEIFVAN
GTQGKLTCKFKSTSTTGGLTSVSWSFQPEGADTTVSFFHYSQGQVYLGNY
PPFKDRISWAGDLDDKDDASINIENMQFIHNGTYICDVKNPPDIVVQPGHI
RLYVVEKENLPVFPVWVVVGIVTAVVLGLTLLISMILAVLYRRKNSKRDY
TGCSTSESLSPVKQAPRKSPSDTEGLVKSLPSGSHQGPVIYAQLDHSGGH
HSDKINKSESVVYADIRKN

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

MPZL1

Full Name

Myelin protein zero like 1

Alternative Names

MPZL1; PZR; PZR_a; PZR_b; PZR1b; MPZL1b; myelin protein zero-like protein 1; immunoglobulin family transmembrane protein; protein zero related; Myelin protein zero like 1

Gene ID

[9019](#)

UniProt ID

[Q95297](#)