

Product Information

MemDX™ Membrane Protein Human TMEM238L (Transmembrane protein 238 like) Full

Length

Cat. No.: **MPC2480K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human TMEM238L membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TMEM238L

Protein Length

Full length

Protein Class

Transporter

TMD

2

Sequence

MLLGSLWGRCHPGRCALFLILALLLDAVGLVLLLLGILAPLSSWDFFIYT
GALILALSLLLWIIWYSLNIEVSPEKLDL

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

TMEM238L

Full Name

Transmembrane protein 238 like

Alternative Names

TMEM238L; FORCP; LINC00675; FOXA1-Regulated Conserved Small Protein; long intergenic non-protein coding RNA 675; putative transmembrane protein LOC100289255; Transmembrane protein 238 like

Gene ID

[100289255](#)

UniProt ID

[A6NJY4](#)