

Product Information

MemDX™ Membrane Protein Human SERINC4 (Serine incorporator 4) Full Length

Cat. No.: **MPC2524K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SERINC4 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SERINC4

Protein Length

Full length

Protein Class

Transporter

TMD

10

Sequence

MVGAKAGPSPGTSLGLAQQHSGGSSVLVKSPFCQVCCCGPAPCASCCHSR
WPSLTASTCSRLFYILLHVGASAI CLLLSRTVVERVWGKTHRIQMPSGL
CAHLFGLSDCPVLSGSGAVYRVCAGTATFHLQAVLLVHLHSPTSPRAQL
HNSFWLLKLLFLLGLCAIAFCIPDEHLFP AWHYIGICGGFAFILLQLVLI
TAFASWNKNWQTGAAQDCSWFLAVLLATLGFYSMAGVGAVLLFHYYTHP
AGCLLNKMLLSLHLCFCGLISFLSIAPCIRLKQPRSGLLQASVISCYIMY
LTF SALSSRPPERVILQGQNHTLCLPGLSKMEPQTPDISLAMLSASIMYA
CVLFACNEASYLAEVFGPLWIVKVYSYEFQKPSLCFCCPETVEADKGQRG
GAARPADQETPPAPPVQVQHLSYNYS AFHFVFFLASLYVMVTLTNWFSYE
GAELEKTFIKGSWATFWVKVASCWACVLLYLGLLLAPLCWPPTQKPQPLI
LRRRRHRIISPDNKYPPV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SERINC4

Full Name

Serine incorporator 4

Alternative Names

SERINC4; Serine incorporator 4

Gene ID

[619189](#)

UniProt ID

[A6NH21](#)