

Product Information

MemDX™ Membrane Protein Human SMIM43 (Small integral membrane protein 43) Full

Length

Cat. No.: **MPC2561K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SMIM43 membrane protein expressed in *E.coli*. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SMIM43

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MEWELNLLLYLALFFFLFLLFLLLFVVIKQLKNSVANTAGALQPGRLSV
HREPWGFSREQAV

Product Description

Expression Systems

E.coli

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

SMIM43

Full Name

Small integral membrane protein 43

Alternative Names

SMIM43; TMEM155; protein TMEM155; transmembrane protein 155; Small integral membrane protein 43

Gene ID

[132332](#)

UniProt ID

[Q4W5P6](#)