

Product Information

MemDX™ Membrane Protein Human YIF1A (Yip1 interacting factor homolog A, membrane trafficking protein) Full Length

Cat. No.: **MPC2669K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human YIF1A membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

YIF1A

Protein Length

Full length

Protein Class

Transporter

TMD

5

Sequence

MAYHSGYGAHGSKHRARAAPDPPLFDDTSGGYSSQPGGYPATGADVAFS
VNHLGDPMANVAMAYGSSIASHGKDMVHKELHRFVSVSKLKYFFAVDTA
YVAKKLGLLVFPYTHQNWEVQYSRDAPLPPRQDLNAPDLYIPTMAFITYV
LLAGMALGIQKRFSPEVLGLCASTALVWVMEVLALLGLYLATVRSCLS
TFHLLAYSGYKYVGMILSVLTGLLFGSDGYVALAWTSSALMYFIVRSR
TAALGPDSMGGPVPRQRLQLYLTGAAAFQPLIYWLTFFHLVR

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

YIF1A

Full Name

Yip1 interacting factor homolog A, membrane trafficking protein

Alternative Names

YIF1A; 54TM; YIF1; YIF1P; FinGER7; protein YIF1A; 54TMp; YIP1-interacting factor homolog A; Yip1p-interacting factor; putative Rab5-interacting protein; putative transmembrane protein 54TMp; Yip1 interacting factor homolog A, membrane trafficking protein

Gene ID

[10897](#)

UniProt ID

[Q95070](#)