

Product Information

MemDX™ Membrane Protein Human ZDHHC13 (Zinc finger DHHC-type palmitoyltransferase

13) Full Length

Cat. No.: **MPC2733K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ZDHHC13 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ZDHHC13

Protein Length

Full length

Protein Class

Transferase

TMD

6

Sequence

MEGPGLGSQCRNHSHPHPPGFGRYGICAHENKELANAREALPLIEDSSN
CDIVKATQYGIFERCKELVEAGYDVRQPDKENVSLLHWAAINNRLDLVKF
YISKGA VVDQLGGDLNSTPLHWAIRQGHLP MVILLQLHGADPTLIDGEGF
SSIH L AVL FQHMP I IAYLISKGQSVNMTDVNGQTPLMLSAHKVIGPEPTG
FLLKFNPSLNVVDKIHQNTPLHWAVAAGNVNAVDKLEAGSSLDIQNVKG
ETPLDMALQNK NQLIIHMLKTEAKMRANQKFRLWRWLQKCELFLLMLSV
ITMWAIGYILDFNSDSWLLKGCLLVTLFFLTSLFPRFLVGYKNLVYLPTA
FLLSSVFWIFMTWIFLFFPDLAGAPFYFSFISIVAFLYFFYKTWATDPG
FTKASEEEKKVNIITLAETGSLDFRTFCTSLIRKPLRSLHCHVCNCCVA
RYDQHCLWTGRCIGFGNHHYIIFLFFLSMVCGWIIYGSFIYLSHCATT
FKEDGLWTYLNQIVACSPWVLYILMLATFHFSWSTFLLNQLFQIAFLGL
TSHERISLQKQSKHMKQTLRLKTPYNLGMQNLADFFQC GCFGLVKPCV
VDWTSQYTMVFHPAREKVLRSV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ZDHHC13

Full Name

Zinc finger DHHC-type palmitoyltransferase 13

Alternative Names

ZDHHC13; HIP14L; HIP3RP; palmitoyltransferase ZDHHC13; DHHC-13; HIP14-related protein; huntingtin interacting protein HIP3RP; huntingtin-interacting protein 14-related protein; probable palmitoyltransferase ZDHHC13; putative MAPK-activating protein PM03; putative NF-kappa-B-activating protein 209; zinc finger DHHC domain-containing protein 13; zinc finger DHHC-type containing 13; zinc finger, DHHC domain containing 13; Zinc finger DHHC-type palmitoyltransferase 13

Gene ID

[54503](#)

UniProt ID

[Q8IUH4](#)