

Product Information

MemDX™ Membrane Protein Human INAFM2 (InaF motif containing 2) Full Length

Cat. No.: **MPC2757K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human INAFM2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

INAFM2

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MKERDAAPAERGKPATYTGDKKAKMAAKTNKKWVRLATVFAYVLSVSLAA
IVLAVYYSLIWQPVGAGTSGGAAGPPPGGSNATGPSGTSGAAAAGPNTTG
SSRREAPRDVPPLQAARPAPPEPPADSPAGPLERPRGPDEDEEETAAAP
GSR

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

INAFM2

Full Name

InaF motif containing 2

Alternative Names

INAFM2; OGU1; LINC00984; putative transmembrane protein INAFM2; inaF-motif-containing protein 2; osteogenesis up-regulated transcript 1; osteogenesis upregulating transcript 1; InaF motif containing 2

Gene ID

[100505573](#)

UniProt ID

[P0DMQ5](#)