

Product Information

MemDX™ Membrane Protein Human C17orf113 (Chromosome 17 open reading frame 113)

Full Length

Cat. No.: **MPC2760K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human C17orf113 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

C17orf113

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MVPPGKKPAGEASNSNKKCKRYFNEHWKEEFTWLDFDYERKLMFCLECRQ
ALVRNKHGKAENAFVTGDNFQRHALLRHVTSGAHRQALAVNQGQPPFEG
QAEGGGACPGlattPASRGVKVELDPAKVAVLTTVYCMakedVPNDRCsa
LLELQRFNLCQALLGTEHGDYYSPRRVRDMQVAIASVLHTEACQRLKASP
YVGLVLDETRDWPESHSLALFATSVSPCDGQPATTFLGSVELQEGEATAG
QLLDILQAFGVSAPKLAWLSSSLPSERLGSVGPQLRATCPLLAELHCLPG
RTDPEPPAYLGQYESILDALFRLHGGPSSHLVPELRAALDLAIDLAGPR
PVPWASLLPVVEAVAEAWPGLVPTLEAAALASPvAGSLALALRQFTFVAF
THLLLDALPSVQKLSVLQAEEDLALLQPLVMAAAASLQAQRGSGGARL
QGFLQELASMDPDASSGRCTYRGVELLGYSEAAVRGLEWLRGSFLDSMRK
GLQDSYPGPSLDAVAFAAIFDPRRYPQAPEELGTHGEGALRVLLRGFAP
AVVRQRALGDFALFKRVVFLGLRGLPRALCTQLACAHSELHELFPDFAAL
AALALALPAGAGLLDKVGRSRELRWVGQSGAGEGRGGMVKIAVDGPPLH
EFDfGLAVEfLESGWGEGFLGSQLT

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

C17orf113

Full Name

Chromosome 17 open reading frame 113

Alternative Names

C17orf113; transmembrane protein C17orf113; Chromosome 17 open reading frame 113

Gene ID

[110806298](#)

UniProt ID

[A0A1B0GUU1](#)