

Product Information

MemDX™ Membrane Protein Human SPTY2D1OS (SPTY2D1 opposite strand) Full Length

Cat. No.: **MPC2761K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SPTY2D1OS membrane protein expressed in *E.coli*. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SPTY2D1OS

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MIVLGWMFFVGLVCYMGTFPELMPPTLKWQERWPVQESKTQLRRRALGED
LLQNHVEGI

Product Description

Expression Systems

E.coli

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw

cycles.

Target

Target Protein

SPTY2D1OS

Full Name

SPTY2D1 opposite strand

Alternative Names

SPTY2D1OS; SPTY2D1-AS1; putative transmembrane protein SPTY2D1OS; SPTY2D1 antisense RNA 1 (non-protein coding); SPTY2D1 opposite strand

Gene ID

[100506540](#)

UniProt ID

[A0A0U1RRN3](#)